

MicroPen®

WHITE LINE




www.micromed.biz




Micromed
Medizintechnik GmbH
Eisenbahnstraße 84
D-78573 Wurmlingen
Tel +49 7461 96 48 55 0
Fax +49 7461 96 48 55 95

SOME SHAPES ARE
SIMPLY PERFECT

SHARK:
PERFECTLY FITTED SHAPE
SINCE 230 MILLION YEARS



MICROPEN®:
MOVING AHEAD - WITH THE
PERFECT FIT



MicroPen®

- High voltage up to 9000 V
- Perfect handling
- Extreme durability
- Modular design

THE MASTER OF ITS LEAGUE

CHOOSE THE BEST!

SEPERABLE MIRCOPEN® AND CABLE

- multifunctional MicroPen® fits on all common electro-surgical units

ELECTRODE-HANDLE MIRCOPEN®

- optimized handling
- aesthetic design

OPTIMUM SYNTHETIC SUBSTANCE

- high electric strength up to 9000 volts
- more than 250 sterilization cycles
- highly resistant to detergent

CONNECTION

- various types of plugs for different manufacturers, RF surgical units
 - 1 Aesculap, Berchtold, Martin
 - 2 ERBE, VIO, ICC und ACC
 - 3 Valleylab 3-PIN-plug
- in various lengths (3m, 5m)

1
order numbers
100-228-340, 3m
100-228-540, 5m

2
order numbers
100-238-340, 3m
100-238-540, 5m

3
order numbers
100-248-340, 3m
100-248-540, 5m

PUSH BUTTON


- lift off stroke can be felt easily
- click function
- high security by durable push button

ELECTRODES

- hexagonal fit for electrodes measuring 4mm and 3/32 inch in diameter
- electrodes can be exchanged easily

LONG DURABILITY

- electronics completely covered with insulating silicone



seperable MicroPen® and cable

THE MASTER OF ITS LEAGUE

- High voltage up to 9000 V
- Perfect handling
- Extreme durability
- Modulary design

MicroPen®

WHITE LINE

12
MONTH
GARANTY*

*MicroPen® comes with a 12 months garanty on purchase